

BEFOUR

MEASURING NEW HEIGHTS

Display Features

Measurements are displayed in bright LED digits in Inches, Feet and Inches, or Centimeters



Inches



Feet and Inches



Centimeters

Flexible Measurement Range

Measures up to 91 inches



DIGITAL HEIGHT ROD: Model HTR-101 Standalone Wall-Mount Digital Stadiometer

Befour's new HTR-101 Height Rod features fast and accurate digital performance at an attractive price. With its user-selectable 51" measurement range, the versatile HTR-101 can measure the height of a small child to that of an NBA player. The smooth glide headpiece features one-button operation, bright LED digits that can be easily seen even at maximum height, and a unique auto-hold function retains the measurement display even as it is raised out of the way.

The HTR-101 stadiometer has an accuracy of ± 0.1 in (0.1 cm) and easily converts between three display formats – inches, feet & inches, and centimeters. With no wires to get in the way, the HTR-101 is ideal for clinics, hospitals, schools, and sport and fitness centers. Powered by (2) "AAA" alkaline batteries (included), and with its automatic turn off, the HTR-101 is good for up to 10,000 measurements.



FEATURES

- Fast and accurate measurement
- View height in inches - feet and inches - centimeters
- Measures up to 91 in (7 ft & 7 in / 230 cm)
- Effortless sliding headpiece
- Flexible measurement range
- Bright LED digits

SPECIFICATIONS

Accuracy	± 0.1 in (0.1 cm)
Measuring Span	51 in (130 cm)
Display	Bright LED Digits
Power	(2) "AAA" batteries (included)
Weight.....	2.5lbs
Rail Dimensions (HxWxD)	57.5 x 2.2 x 0.5 in (146 x 5.6 x 1.3 cm)
Headpiece	Extends 8 inches from wall (fixed)
Construction.....	Anodized aluminum rail
Warranty	3 years

BEFOUR, INC

102 Progress Drive
Saukville, WI 53080
Phone: 800-367-7109
Fax: 262-284-5966
Email: sales@befour.com
www.befour.com

