

BEFOUR

PRECISION SCALES

HEALTH & FITNESS SCALE

Quick, reliable, and extremely accurate in high-volume environments...

Model FS-0900 Free Standing Scale

Stand-on scales are ideal for locker rooms, fitness clubs, medical facilities, and personal use.



With its plug-in design, the FS-0900 is a general purpose scale that provides quick, accurate, and reliable weight measurements. Designed with functional basics in mind, our stand-on scale is ideal for locker rooms, fitness facilities, and personal use.

- AC plug-in style
- Large platform for greater stability
- Pin point accuracy of ± 0.1 lb.
- Huge 500 pound capacity
- Rugged, lightweight, all metal construction for years of service
- Easy assembly

MOTIONLOCK

MODEL FS-0900

Scale Height.....	42"
Platform Dimensions	16" x 18" x 1 $\frac{3}{4}$ "
Weight Capacity	500 lbs. or 225 kg.
Accuracy / Resolution.....	± 0.1 lb. (± 0.1 kg.)
Power Source.....	120V AC
Readout.....	Bright $\frac{3}{4}$ " LED digits
System Weight	17 lbs.
Warranty	3 Years

FEATURES/SPECIFICATIONS:

Push Button Calibration: Set calibration via keypad entry without the need of tools

Auto Zero: Scale automatically returns to zero after each weighing

Auto Off, Tare

Unique MotionLock Technology: Incredibly accurate measurements even for unsteady users

Weight Measurements: In pounds or kilograms



Proudly Made in the USA

LED DISPLAY FEATURES:



- Pinpoint Accuracy: ± 0.1 lb. makes it extremely accurate when weight is a factor
- LED Display: Bright $\frac{3}{4}$ " LED digits can be easily seen
- Power source: AC plug-in style
- Push Button Calibration: Set calibration via keypad entry; no tools needed
- MotionLock Technology: Accurate measurements even for the most unsteady user
- Auto-Zero: The scale automatically returns to zero after each weighing
- Auto off, tare functions
- Weight Measurement: In pounds or kilograms



Befour Inc.
102 Progress Drive
Saukville, WI 53080

Phone
800-367-7109

Web
www.befour.com

Fax
262-284-5966

E-mail
sales@befour.com