

# BEFOUR

## PRECISION SCALES



Proudly Made in the USA

## HEALTH & FITNESS SCALE

Quick, reliable, and extremely accurate in high-volume environments...

### Model FS-0961 Pro-BMI Stand-on Scale

Stand-on scales are ideal for locker rooms, fitness clubs, medical facilities, and personal use.



Battery power allows user to place scale anywhere-no outlet or re-charging needed!  
Efficient design allows batteries to last over 80,000 weighings.  
Quick 2 second weight measurement, plus one button weight recall and reweigh in either pounds or kilograms

- Quickly displays Body Mass Index
- Pin point accuracy of  $\pm 0.1$  lb.
- Large platform for greater stability
- Programmable to your specific needs
- Rugged, lightweight, all-metal construction
- Easy assembly

#### MODEL FS-0961

|                            |                                |
|----------------------------|--------------------------------|
| Scale Height.....          | 42"                            |
| Platform Dimensions.....   | 16" x 18" x 1 $\frac{3}{4}$ "  |
| Weight Capacity.....       | 500 lbs. or 225 kg.            |
| Accuracy / Resolution..... | $\pm 0.1$ lb. ( $\pm 0.1$ kg.) |
| Power Source.....          | 6 C batteries                  |
| Readout.....               | 1" LCD digits                  |
| System Weight.....         | 17 lbs.                        |
| Warranty.....              | 3 Years                        |



### FEATURES/SPECIFICATIONS:

**Body Mass Index (BMI) Calculator:** Dual-mode BMI calculator allows user - or assistant-operation and provides a quick and easy BMI calculation.

**Unique MotionLock Technology:** Incredibly accurate weight measurements even for unsteady users.

**Push-button Calibration:** Set calibration via keypad entry with no tools needed.

**"Always-on" Option:** For unequalled speed of measurement, set up the console so that no buttons need to be pressed for normal weighings.

**Programmable setup:** Many functions are programmable to setup to your specific needs.

### LCD DISPLAY FEATURES:



- Alkaline batteries provide over 80,000 weighings
- BMI calculator with programmable clothing allowance
- "Always-On" option
- Weight reads in pounds or kilograms
- Programmable pound/kilogram lockout
- Programmable auto-off:
- 12-seconds to 25-minutes
- Auto zero
- Weight recall and reweigh
- Audible sound indicates weight is "locked-in"
- Durable anodized aluminum construction
- Large 1" easy-to read digits



Befour Inc.

102 Progress Drive  
Saukville, WI 53080

Phone

800-367-7109

Web

www.befour.com

Fax

262-284-5966

E-mail

sales@befour.com