

HEALTH & FITNESS SCALE



Proudly Made in the USA

Quick, reliable, and extremely accurate in high-volume environments...

Model FS-0961 Pro-BMI Stand-on Scale

Stand-on scales are ideal for locker rooms, fitness clubs, medical facilities, and personal use.



Battery power allows user to place scale anywhere-no outlet or re-charging needed!
Efficient design allows batteries to last over 80,000 weighings. Quick 2 second weight measurement, plus one button weight recall and reweigh in either pounds or kilograms

- Quickly displays Body Mass Index
- Pin point accuracy of ± 0.1 lb.
- Large platform for greater stability
- Programmable to your specific needs
- Rugged, lightweight, all-metal construction
- Easy assembly

MODEL FS-0961

Scale Height.....	42"
Platform Dimensions.....	16" x 18" x 1 $\frac{3}{4}$ "
Weight Capacity.....	500 lbs. or 225 kg.
Accuracy / Resolution.....	± 0.1 lb. (± 0.1 kg.)
Power Source.....	6 C batteries
Readout.....	1" LCD digits
System Weight.....	17 lbs.
Warranty.....	3 Years



FEATURES/SPECIFICATIONS:

Body Mass Index (BMI) Calculator: Dual-mode BMI calculator allows user - or assistant-operation and provides a quick and easy BMI calculation.

Unique MotionLock Technology: Incredibly accurate weight measurements even for unsteady users.

Push-button Calibration: Set calibration via keypad entry with no tools needed.

"Always-on" Option: For unequalled speed of measurement, set up the console so that no buttons need to be pressed for normal weighings.

Programmable setup: Many functions are programmable to setup to your specific needs.

LCD DISPLAY FEATURES:



- Alkaline batteries provide over 80,000 weighings
- BMI calculator with programmable clothing allowance
- "Always-On" option
- Weight reads in pounds or kilograms
- Programmable pound/kilogram lockout
- Programmable auto-off:
- 12-seconds to 25-minutes
- Auto zero
- Weight recall and reweigh
- Audible sound indicates weight is "locked-in"
- Durable anodized aluminum construction
- Large 1" easy-to read digits

