



## HEALTHCARE SCALE COMPARISON GUIDE

...the important information you need at a glance

										
	<b>MX115</b> 550 lbs.	<b>MX160</b> 750 lbs.	<b>MX170</b> 1000 lbs.	<b>MX305</b> 750 lbs.	<b>MX306</b> 750 lbs.	<b>MX308</b> 750 lbs.	<b>MX310</b> 1000 lbs.	<b>MX805</b> 750 lbs.	<b>MX808</b> 750 lbs.	<b>MX810</b> 1000 lbs.
Weight Capacity	550 lbs.	750 lbs.	1000 lbs.	750 lbs.	750 lbs.	750 lbs.	1000 lbs.	750 lbs.	750 lbs.	1000 lbs.
Platform Size (width x depth)	13" x 13"	16" x 18"	24" x 20"	16" x 18"	16" x 18"	16" x 18"	22" x 25"	16" x 18"	16" x 18"	22" x 25"
Scale/System Weight	12 pounds	14 pounds	25 pounds	17 pounds	31 pounds	31 pounds	46 pounds	17 pounds	33 pounds	46 pounds
Resolution	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound
Accuracy	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound
Portability	Carry	Carry	Carry Handle		Tilt & Roll	Tilt & Roll	Tilt & Roll		Tilt & Roll	Tilt & Roll
Handrail					Standard	Standard	Standard		Standard	Standard
Integrated Height Rod								◆	◆	◆
Alkaline Battery Power	◆	◆	◆	◆	◆	◆	◆	◆	◆	◆
Body Mass Index (BMI)	◆	◆	◆	◆	◆	◆	◆	◆	◆	◆
Stand-on	◆	◆	◆	◆	◆	◆	◆	◆	◆	◆
Chair			◆							

							
	<b>MX202</b> 45 lbs.	<b>MX420</b> 750 lbs.	<b>MX450</b> 1000 lbs.	<b>MX450D</b> 1000 lbs.	<b>MX480</b> 1000 lbs.	<b>MX480D</b> 1000 lbs.	<b>MX490D</b> 1000 lbs.
Weight Capacity	45 lbs.	750 lbs.	1000 lbs.	1000 lbs.	1000 lbs.	1000 lbs.	1000 lbs.
Platform Size (width x depth)	26" x 16 1/2"	Variable	28" x 32"	28" x 32"	32" x 34"	32" x 34"	32" x 41"
Ramps		Integral	Single	Dual	Single	Dual	Dual
Scale/System Weight	27 pounds	24 pounds	65 pounds	65 pounds	87 pounds	87 pounds	100 pounds
Resolution	1 gram	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound
Accuracy	± 1 gram	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound
Portability	Handles	Carry	Tilt & Roll	Tilt & Roll	Fold & Roll	Fold & Roll	Fold & Roll
Handrail			Standard	Standard	Standard	Standard	Standard
Alkaline Battery Power	◆	◆	◆	◆	◆	◆	◆
Body Mass Index (BMI)		◆	◆	◆	◆	◆	◆
Stand-on			◆	◆	◆	◆	◆
Chair		◆	◆	◆	◆	◆	◆
Wheelchair		◆	◆	◆	◆	◆	◆

**Befour's precision Medical scales have  Technology, the medical console displays have 1" LCD digits and ALL come with a \*\*3 year Warranty\*\***