



HEALTHCARE SCALE COMPARISON GUIDE

...the important information you need at a glance

										
	MX115 550 lbs.	MX160 750 lbs.	MX170 1000 lbs.	MX305 750 lbs.	MX306 750 lbs.	MX308 750 lbs.	MX310 1000 lbs.	MX805 750 lbs.	MX808 750 lbs.	MX810 1000 lbs.
Weight Capacity	550 lbs.	750 lbs.	1000 lbs.	750 lbs.	750 lbs.	750 lbs.	1000 lbs.	750 lbs.	750 lbs.	1000 lbs.
Platform Size (width x depth)	13" x 13"	16" x 18"	24" x 20"	16" x 18"	16" x 18"	16" x 18"	22" x 25"	16" x 18"	16" x 18"	22" x 25"
Scale/System Weight	12 pounds	14 pounds	25 pounds	17 pounds	31 pounds	31 pounds	46 pounds	17 pounds	33 pounds	46 pounds
Resolution	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound
Accuracy	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound
Portability	Carry	Carry	Carry Handle		Tilt & Roll	Tilt & Roll	Tilt & Roll		Tilt & Roll	Tilt & Roll
Handrail					Standard	Standard	Standard		Standard	Standard
Integrated Height Rod								◆	◆	◆
Alkaline Battery Power	◆	◆	◆	◆	◆	◆	◆	◆	◆	◆
Body Mass Index (BMI)	◆	◆	◆	◆	◆	◆	◆	◆	◆	◆
Stand-on	◆	◆	◆	◆	◆	◆	◆	◆	◆	◆
Chair			◆							

										
	MX202 45 lbs.	MX308CHR 500 lbs.	MX420 750 lbs.	MX450 1000 lbs.	MX450D 1000 lbs.	MX480 1000 lbs.	MX480D 1000 lbs.	MX490D 1000 lbs.	VS0900 600 lbs.	
Weight Capacity	45 lbs.	500 lbs.	750 lbs.	1000 lbs.	1000 lbs.	1000 lbs.	1000 lbs.	1000 lbs.	600 lbs.	
Platform Size (width x depth)	26" x 16 1/2"	21" x 17" seat	Variable	28" x 32"	28" x 32"	32" x 34"	32" x 34"	32" x 41"	42" x 22"	
Ramps			Integral	Single	Dual	Single	Dual	Dual		
Scale/System Weight	27 pounds	64 pounds	24 pounds	65 pounds	65 pounds	87 pounds	87 pounds	100 pounds	46 pounds	
Resolution	1 gram	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	0.1 pound	
Accuracy	± 1 gram	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	
Portability	Handles	Tilt & Roll	Carry	Tilt & Roll	Tilt & Roll	Fold & Roll	Fold & Roll	Fold & Roll		
Handrail		Standard		Standard	Standard	Standard	Standard	Standard		
Alkaline Battery Power	◆	◆	◆	◆	◆	◆	◆	◆	◆	
Body Mass Index (BMI)		◆	◆	◆	◆	◆	◆	◆	◆	
Stand-on		Converts to Stand On		◆	◆	◆	◆	◆	◆	
Chair		◆	◆	◆	◆	◆	◆	◆	◆	
Wheelchair			◆	◆	◆	◆	◆	◆		

Befour's precision Medical scales have **MOTIONLOCK** Technology, the Medical Console

Displays have 1" LCD digits and ALL come with a ****3 year Warranty****

Proudly Made in the USA!!!