

Professional Healthcare Scales

BEFOUR

1-800-367-7109

www.befour.com

CX7 CONSOLE FEATURES:



- Displays weight, height, and BMI on a single screen
- EMR connectivity via serial output or optional wireless module
- ECO Weigh battery system provides 100,000 weighings per set of (4) D-cell alkaline batteries (included)
- 3-segment battery life indicator
- Large 1" LCD digits
- No AC outlet or recharging needed



Befour products are made with pride in the USA!


Befour, Inc.

102 Progress Drive
Saukville, WI 53080
Toll free: 800-367-7109
Phone: 262-284-5150
Fax: 262-284-5966

Model **MX115**

Portable Platform Scale with BMI

This is a truly versatile scale ideal where precise weight measurements are required. It's quality and ease of use are perfect for visiting nurses, weight management clinics, and any situation where high capacity and lightweight are essential. Convenient alkaline battery power allows user to place scale anywhere – no outlet or re-charging needed!



FEATURES/SPECIFICATIONS

High capacity: Huge 550 pound capacity.

Maximum Portability: All-aluminum design weighs only 12 pounds!

Extreme Accuracy: ± 0.1 lbs.

Body Mass Index: BMI calculator

MotionLock Technology: Accurate weight even with an unsteady user.

LB/KG Lockout: Quickly change setup to disable LB/KG button so scale will only weigh in selected unit of measure.

New CX7 Console: Displays only the information you need - weight only or weight, height, and BMI.

Versatile Console: Mounting brackets (attached) allow positioning on wall or angled on tabletop.

Platform Dimensions	13" x 13" x 1.75"
System Weight	12 lbs.
Weight Capacity	550 lbs. / 250 kg
Resolution	0.1 lb. / (0.05 kg.)
Accuracy	± 0.1 lb. / (0.05 kg.)
Console Mounting	Wall, table brackets (included)
Power Source	4 D-cell alkaline batteries
Connectivity	USB Serial or optional wireless
Warranty	3 Years