

# Professional Healthcare Scales

**BEFOUR**

**1-800-367-7109**

**www.befour.com**

## CX7 CONSOLE FEATURES:



- Displays weight, height, and BMI on a single screen
- EMR connectivity via serial output or optional wireless module
- ECO Weigh battery system provides 100,000 weighings per set of (4) D-cell alkaline batteries (included)
- 3-segment battery life indicator
- Large 1" LCD digits
- No AC outlet or recharging needed



**Befour products are made with pride in the USA!**

### Befour, Inc.

102 Progress Drive  
Saukville, WI 53080  
Toll free: 800-367-7109  
Phone: 262-284-5150  
Fax: 262-284-5966

## Model MX170

### Portable Bariatric Scale with BMI

Ultra-rugged design with no detectible flex underfoot!  
Oversized, rock-solid platform incorporates triple-walled extruded aluminum framework, full-perimeter tensile reinforcement. Alkaline battery allows use anywhere.

Oversized platform and high capacity offer an ideal scale for bariatric clinics.



### FEATURES/SPECIFICATIONS

**High capacity:** Huge 1000 pound capacity.

**Oversize Platform:** For user safety and security.

**Easy transport:** Lightweight all-metal construction with built-in handle.

**Versatile Console:** Console stores underneath platform for portability. Mounting brackets (attached) allow positioning on wall or angled on tabletop.

**MotionLock Technology:** Accurate weight even with an unsteady user.

**LB/KG Lockout:** Quickly change setup to disable LB/KG button so scale will only weigh in selected unit of measure.



Built-in handle and recessed storage cavity securely hold readout during storage and extended transport.

Platform Dimensions	24" x 20" x 2"
System Weight	25 lbs.
Weight Capacity	1000 lbs. / 450 kg
Resolution	0.1 lb. / (0.05 kg.)
Accuracy	± 0.1 lb. / (0.05 kg.)
Console Mounting	Wall, table brackets (included)
Power Source	4 D-cell alkaline batteries
Connectivity	USB Serial or optional wireless
Warranty	3 Years

**BEFOUR**