

# Professional Healthcare Scales

**BEFOUR**

**1-800-367-7109**  
**www.befour.com**

**Model CRT202**  
**OPTIONAL CART**  
**FOR NEONATAL SCALE**



MX202 scale becomes a self-contained mobile weighing station when used with optional cart.

- Two 24" w x 18" d shelves with 1/4" safety lip.
- Lockable cabinet with built-in handle
- Four 4" casters (two locking)
- 200 pound capacity
- Rugged thermoplastic cart, painted steel cabinet
- Easy 10 minute assembly



*Befour scales are made with pride in the USA!*

**Befour, Inc.**  
 102 Progress Drive  
 Saukville, WI 53080  
 Toll free: 800-367-7109  
 Phone: 262-284-5150  
 Fax: 262-284-5966

## Model MX202

### Neonatal/Pediatric Scale with 1 Gram resolution

Designed for daily and serial weighing of neonates and infants. Heavy-duty all-metal portable unit with a large removable acrylic cradle with length rule. A standard 1 gram accuracy and large 45 pound (20 kg.) capacity means this single scale can easily handle all infant care weighing needs. The accurate and repeatable readings are ideal for determining breast-feeding intake or urine output by diaper weighing.

Engineered for maximum mobility and safety, the scale has two large carrying handles. Six standard C-cell alkaline batteries (included) provide over 80,000 measurements with no outlet or recharging needed.



#### FEATURES/SPECIFICATIONS

**1 gram resolution:** Standard with this scale for extremely accurate weighing.

**MotionLock Technology:** Accurately calculates the weight of active babies.


**No set-up required:** Arrives fully assembled and ready to use.

**LB/GM Lockout:** Quickly change setup to disable LB/GM button so scale will only weigh in selected unit of measure: gm, kg, or lb/oz.




**Optional 4-sided acrylic cradle**

|                            |  |
|----------------------------|--|
| Dimension with Handles ... | 26" x 16 1/2" x 7"   |
| Acrylic Cradle Size .....  | 24" x 15" x 3"   |
| System Weight .....        | 27 lbs   |
| Weight Capacity .....      | 20,000 grams / 45 lbs.   |
| Resolution .....           | 1 gram (to 10,000 grams)<br>5 grams (10,001 to 20,000 grams)<br>0.1 oz over entire range       |
| Accuracy .....             | ± 1 gram (to 10,000 grams)<br>± 5 grams (10,001 to 20,000 grams)<br>± 0.1 oz over entire range |
| Power Source .....         | 6 C-cell alkaline batteries  |
| Display .....              | 1" LCD   |
| Warranty .....             | 3 years  |



Grams



Kilograms



Pounds/Ounces

Three display options: grams, kilograms, or pounds/ounces.



1 gram resolution for diaper weighing!