

Professional Healthcare Scales

BEFOUR

1-800-367-7109

www.befour.com

CX7 CONSOLE FEATURES:



- Displays weight, height, and BMI on a single screen
- EMR connectivity via serial output or optional wireless module
- ECO Weigh battery system provides 100,000 weighings per set of (4) D-cell alkaline batteries (included)
- 3-segment battery life indicator
- Large 1" LCD digits
- No AC outlet or recharging needed




Befour products are made with pride in the USA!

Befour, Inc.

102 Progress Drive
Saukville, WI 53080
Toll free: 800-367-7109
Phone: 262-284-5150
Fax: 262-284-5966

Model **MX305** Column Scale with BMI

Attractive light weight, all-aluminum design features a rugged cast platform and anodized column and console housing for years of dependable operation. Easy viewing 42" display height and Body Mass Index (BMI) calculator.



FEATURES/SPECIFICATIONS

High capacity: Now with huge 750 pound capacity.

Large Stable Platform: Recessed non-skid mat and low step-up height.

Built to last: Rugged, non-corroding all-aluminum design weighs only 17 lbs.

MotionLock Technology: Accurate weight even with an unsteady user.

LB/KG Lockout: Quickly change setup to disable LB/KG button so scale will only weigh in selected unit of measure.

New CX7 Console: Displays only the information you need - weight only or weight, height, and BMI.



Scale Height	42"
Platform Dimensions	16" x 18" x 1.75"
System Weight	17 lbs.
Weight Capacity	750 lbs. / 340 kg
Resolution	0.1 lb. / (0.05 kg.)
Accuracy	± 0.1 lb. / (0.05 kg)
Power Source	4 D-cell alkaline batteries (included)
Connectivity	USB Serial output or optional wireless
Warranty	3 Years