

Professional Healthcare Scales

BEFOUR

1-800-367-7109
www.befour.com

CX7 CONSOLE FEATURES:



- Displays weight, height, and BMI on a single screen
- EMR connectivity via serial output or optional wireless module
- ECO Weigh battery system provides 100,000 weighings per set of (4) D-cell alkaline batteries (included)
- 3-segment battery life indicator
- Large 1" LCD digits
- No AC outlet or recharging needed



Sturdy handgrip for patient stability once on the scale.



Befour scales are made with pride in the USA!

Befour, Inc.
102 Progress Drive
Saukville, WI 53080
Toll free: 800-367-7109
Phone: 262-284-5150
Fax: 262-284-5966

Model MX306

Tilt & Roll Space-Saving Handgrip Scale

The MX306 scale is designed to provide the security of a patient handgrip, but take up minimal floor space. Weight is not affected when patient grabs the bar. The scale is easy to move, just tilt back and roll to the next location. Ideal for hospitals, doctor's offices, and extended care facilities. Convenient battery power means the scale can be placed anywhere.

**Minimal space -
Maximum security**

FEATURES/SPECIFICATIONS

- Stable Handgrip:** Provides security without taking up space.
- Maximum mobility:** Heavy duty wheels and lightweight design make for easy transport.
- Built to last:** Rugged all-metal construction.
- MotionLock Technology:** Accurate weight even with an unsteady user.
- LB/KG Lockout:** Quickly change setup to disable LB/KG button so scale will only weigh in selected unit of measure.
- CX7 console:** Displays only the information you need - weight only or weight, height, and BMI.
- High capacity:** 750 lb. capacity meets most needs!



Easy transport: tilt and roll on heavy-duty 3" wheels.



**EMR
Ready!**

Requires only
17" x 22" of
floor space.

Platform Dimensions	16" x 18" x 1.75"
Scale Height	42"
System Weight	31 lbs.
Weight Capacity	750 lbs. / 340 kg
Resolution	0.1 lb. / (0.05 kg.)
Accuracy	± 0.1 lb. / (0.05 kg.)
Power Source	4 D-cell alkaline batteries
Connectivity	USB Serial or optional wireless
Warranty	3 Years

BEFOUR