

Professional Healthcare Scales

BEFOUR

1-800-367-7109

www.befour.com

CX7 CONSOLE FEATURES:



- Displays weight, height, and BMI on a single screen
- EMR connectivity via serial output or optional wireless module
- ECO Weigh battery system provides 100,000 weighings per set of (4) D-cell alkaline batteries (included)
- 3-segment battery life indicator
- Large 1" LCD digits
- No AC outlet or recharging needed



Befour scales are made with pride in the USA!

Befour, Inc.

102 Progress Drive
Saukville, WI 53080
Toll free: 800-367-7109
Phone: 262-284-5150
Fax: 262-284-5966

Model MX308CHR

Convertible Chair Scale


Unique convertible design gives you two scales in one: a wide body chair scale and a handrail scale. Conversion takes only seconds and requires no tools. First-class seating with 3" padded seat bottom and 25" span between armrests. Easy mobility: simply tilt back and roll on heavy-duty wheels.

FEATURES/SPECIFICATIONS

- Extra Wide Chair:** Accommodates larger individuals.
- Anti-Microbial Vinyl Covering:** For easy cleaning.
- Fast Conversions:** Requires seconds to convert to handrail scale mode.
- Slide-Out Footrest:** For patient comfort.
- MotionLock Technology:** Accurate weight even with an unsteady user.
- LB/KG Lockout:** Quickly change setup to disable LB/KG button so scale will only weigh in selected unit of measure.



EMR Ready!



Slide-out footrest.



MX308CHR in handrail scale mode. (Converts in seconds.)



Chair Scale Dimensions	29"w x 27"d x 42"h
Span Between Armrests	25"
Wheels	3" wheels
System Weight	64 lbs.
Weight Capacity	500 lbs. / 220 kg
Resolution	0.1 lb. / (0.05 kg.)
Accuracy	± 0.1 lb. / (0.05 kg.)
Power Source	4 D-cell alkaline batteries
Connectivity	USB Serial or optional wireless
Warranty	3 Years

BEFOUR