

Professional Healthcare Scales

BEFOUR

1-800-367-7109

www.befour.com

CX7 CONSOLE FEATURES:



- Displays weight, height, and BMI on a single screen
- EMR connectivity via serial output or optional wireless module
- ECO Weigh battery system provides 100,000 weighings per set of (4) D-cell alkaline batteries (included)
- 3-segment battery life indicator
- Large 1" LCD digits
- No AC outlet or recharging needed



Befour scales are made with pride in the USA!

Befour, Inc.

102 Progress Drive
Saukville, WI 53080
Toll free: 800-367-7109
Phone: 262-284-5150
Fax: 262-284-5966

Model MX308

Tilt & Roll Handrail Scale with BMI

Designed to maximize your health care investment, this scale provides the security of a wrap-around handrail with ease of tilt & roll mobility. Ideal for hospitals, doctor's offices, and extended care facilities. Convenient battery power means the scale can be placed anywhere with no cord to trip over.



FEATURES/SPECIFICATIONS

Super-wide wrap-around handrail: Provides maximum security with full 25" between handgrips.

Maximum mobility: Heavy duty wheels and lightweight design make for easy transport.

Built to last: Rugged all-metal construction.

MotionLock Technology: Accurate weight even with an unsteady user.

LB/KG Lockout: Quickly change setup to disable LB/KG button so scale will only weigh in selected unit of measure.

New CX7 console: Displays only the information you need - weight only or weight, height, and BMI.

High capacity: Now with 750 lb. capacity!



Easy transport: tilt and roll on heavy-duty 3" wheels.

Platform Dimensions16" x 18" x 1.75"
Span Between Handgrips 25"
Scale Height 42"
System Weight 31 lbs.
Weight Capacity 750 lbs. / 340 kg
Resolution 0.1 lb. / (0.05 kg.)
Accuracy ± 0.1 lb. / (0.05 kg.)
Power Source 4 D-cell alkaline batteries
Connectivity USB Serial or optional wireless
Warranty 3 Years

BEFOUR