

# Professional Healthcare Scales

**BEFOUR**

**1-800-367-7109**

**www.befour.com**

## CX7 CONSOLE FEATURES:



- Displays weight, height, and BMI on a single screen
- EMR connectivity via serial output or optional wireless module
- ECO Weigh battery system provides 100,000 weighings per set of (4) D-cell alkaline batteries (included)
- 3-segment battery life indicator
- Large 1" LCD digits
- No AC outlet or recharging needed



**Befour scales are made with pride in the USA!**

## Befour, Inc.

102 Progress Drive  
Saukville, WI 53080  
Toll free: 800-367-7109  
Phone: 262-284-5150  
Fax: 262-284-5966

## Model MX420

# Adjustable Width Portable Wheelchair Scale

An excellent choice for the facility that occasionally needs to weigh people in wheelchairs. The dual rails adjust in width to accommodate most wheelchairs. For mobile patients, the scale makes a convenient chair scale - simply place a chair on the rails and you are ready to weigh. Also great for the visiting nurse - it is lightweight and easily fits in the trunk of a car.

When folded for transport, the console and scale store inside carry case (included) and the entire package can be stored in a corner of a closet.

Note: The MX420 is recommended for use on hard flooring surfaces only.



**Width adjusts from 16-3/4" to 35-1/2"**

## FEATURES/SPECIFICATIONS

**Adjustable-width rails:** Conveniently slide to accommodate a variety of wheelchair widths.

**Ergonomic folding design:** Folds in half for easy transport..

**Battery powered:** Allows maximum portability.

**MotionLock Technology:** Accurate weight even with an unsteady user.

**LB/KG Lockout:** Quickly change setup to disable LB/KG button so scale will only weigh in selected unit of measure.

**New Smart View Display:** Displays only the information you need - weight only or weight, height, and BMI.



Scale folds up and stores in nylon carry case for easier transport and storage.

Adjustable Dimensions .....	31" x 16 3/4" adjusts to 35 1/2"
Rail Width .....	3 1/2"
System Weight .....	24 lbs.
Weight Capacity .....	750 lbs. / 340 kg
Resolution .....	0.1 lb. / (0.05 kg.)
Accuracy .....	± 0.1 lb. / (0.05 kg.)
Power Source .....	4 D-cell alkaline batteries
Connectivity .....	USB Serial or optional wireless
Warranty .....	3 Years