

# Professional Healthcare Scales

**BEFOUR**

**1-800-367-7109**

**www.befour.com**

## CX7 CONSOLE FEATURES:



- Displays weight, height, and BMI on a single screen
- EMR connectivity via USB serial output or optional wireless module
- ECO Weigh battery system provides 100,000 weighings per set of (4) D-cell alkaline batteries (included)
- 3-segment battery life indicator
- Large 1" LCD digits
- No AC outlet or recharging needed



**Befour products are made with pride in the USA!**

**Please Note:** The MX700 scale is a custom order product that will require the customer to enlist the aid of a contractor for proper scale pit installation.

## Model MX700


### In-Floor Wheelchair Scale with Handrail



**EMR Ready!**

The oversized platform is flush with the floor, eliminating the need for a ramp. For patients able to stand, the integrated handrail can be grasped and still provide an accurate measurement. MotionLock Technology automatically calculates the weight of even the most unsteady patient.

The standard model has the console mounted on the handrail and no external or floor wiring needed. Can also be custom-ordered with a wall mounted console.



## FEATURES/SPECIFICATIONS

**Flush with flooring:** No need for a ramp.

**Self-contained design:** Platform, handrail, and display are integrated to lift out of the pit as a single unit with no wire routing through the floor or wall.

**Standard Integrated handrail:** Patient's weight is accurately calculated even when using the handrail.

**Generous platform size:** 32" x 50" platform accommodates all wheelchairs.

**Pivoting console:** Pivots 320° for optimal viewing.

**MotionLock Technology:** Accurate weight even with an unsteady user.

**Stainless steel:** Platform cover and pit frame.

**Detailed pit instructions:** To aid contractors.



Heavy-duty handrail is standard with the MX700. The display console rotates for optimal view. Unique self-contained design requires no special wire routing.

### Model MX700

Scale Height .....	42"
Platform Dimensions .....	32" x 50"
System Weight .....	120 lbs.
Weight Capacity .....	1000 lbs. / 450 kg
Resolution .....	0.1 lb. / (0.05 kg.)
Accuracy .....	± 0.1 lb. < 550 lbs > ± 0.1%
Power Source .....	4 D-cell alkaline batteries
Connectivity .....	USB Serial or optional wireless
Warranty .....	3 Years