

Professional Healthcare Scales

BEFOUR

1-800-367-7109

www.befour.com

CX7 CONSOLE FEATURES:



- Displays weight, height, and BMI on a single screen
- EMR connectivity via serial output or optional wireless module
- ECO Weigh battery system provides 100,000 weighings per set of (4) D-cell alkaline batteries (included)
- 3-segment battery life indicator
- Large 1" LCD digits
- No AC outlet or recharging needed



Integrated digital height rod quickly and accurately measures height. Height rod rotates 360° for easier positioning of head piece.

Befour, Inc.

102 Progress Drive
Saukville, WI 53080
Toll free: 800-367-7109
Phone: 262-284-5150
Fax: 262-284-5966

Contract No. V797D-30215
Certified Small Business

Model MX805

Measurement Station Column Scale

The MX805 measuring station offers an attractive light weight, all-aluminum design featuring a rugged cast platform and anodized column and console housing for years of dependable operation. Easy viewing 42" display height.

With its' integrated digital height rod, this scale quickly determines weight, height, and BMI.

FEATURES/SPECIFICATIONS

- High capacity:** Now with huge 750 pound capacity.
- Large Stable Platform:** Recessed non-skid mat and low step-up height.
- Digital Height Measurement:** Simply press HEIGHT and raise the foldable headpiece. Height, in selectable units, displays automatically.
- MotionLock Technology:** Accurate weight even with an unsteady user.
- Built to last:** Rugged, non-corroding all-aluminum design weighs only 17 lbs.
- LB/KG Lockout:** Quickly change setup to disable LB/KG button so scale will only weigh in selected unit of measure.
- New CX7 Console:** Displays only the information you need - weight only or weight, height, and BMI.

EMR Ready!



Scale Height	42"
Platform Dimensions	16" x 18" x 1.75"
System Weight	17 lbs.
Weight Capacity	750 lbs. / 340 kg
Resolution	0.1 lb. / (0.05 kg.)
Accuracy	± 0.1 lb. / (0.05 kg.)
Power Source	4 D-cell alkaline batteries
Connectivity.....	USB Serial or optional wireless
Height Rod Range.....	38 - 88" / (97 - 225 cm)
Height Rod Accuracy	± 0.1 in. / (0.2 cm.)
Height Display Formats	In, Ft & In, cm
Warranty	3 Years



Befour scales are made with pride in the USA!