

BEFOUR

PRECISION SCALES



Provides incredibly accurate weight measurements for unsteady users.

DISPLAY FEATURES: PS-7700



- Easy to use Body Mass Index (BMI) calculation with programmable clothing allowance
- Pinpoint accuracy of ± 0.1 lb.
- The console is powered by six (6) standard C-cell alkaline batteries (included)
- Large 1" LCD digits
- Auto-off time from 12 seconds to 25 minutes, audible tone when weight locks in, lb/kg lockout mode
- Weigh in pounds or kilograms
- 3-segment battery life indicator



Proudly Made in the USA

Constructed with oversized platforms and triple-reinforced undercarriages...

BARIATRIC SCALE

High-capacity, portable Pro-BMI scales are designed with precision and heavy weighing in mind.

Model PS-7700 Pro-BMI

With its generous 24" x 20" platform and a 900 lb. weight capacity, the PS-7700 is ideal for bariatric professionals, home health care, weight-lifting competitions, and athletic weigh-ins. The low-profile platform sits just 2" off the ground, making it easy to step on for individuals of all sizes.



Exclusive "Plug & Weigh" replaceable interconnect cable

- **Weight Capacity:** 900 lb.
- **Oversized Platform:** 20" x 24" x 2"
- **Battery Powered:** Operates on six (6) C-cell alkaline batteries (included)
- **Body Mass Index (BMI):** Quick and easy BMI calculation
- **Accuracy:** ± 0.1 lb. / 0.1 kg
- **Console Display:** 1" LCD digits
- **3 year warranty**

MODEL PS-7700

Platform Dimensions.....	20" x 24" x 2"
Weight Capacity.....	900 lbs. or 405 kg.
Accuracy / Resolution.....	± 0.1 lb. (± 0.1 kg.)
Power Source.....	.6 C batteries
Readout.....	1" LCD digits
System Weight.....	.25 lbs.
Warranty.....	3 Years

FEATURES/SPECIFICATIONS

Body Mass Index (BMI) Calculator: Calculates body mass index in seconds in two easy-to-use modes. Standard mode, used in most situations, simply follow the instructions on the console to input height and display BMI. Assistant mode allows a fitness trainer or health professional to enter a person's height before weighing, eliminating any need to access the keypad while the client is standing on the scale

Unique MotionLock Technology: Incredibly accurate measurements even for unsteady users

High Capacity Tare: Resets scale to zero, even with weight on platform (i.e. mother-child function)

Push Button Calibration: Set calibration via keypad entry without the need of tools

Programmable Set Up: Tailor the console for your specific needs including auto-off time, audible beep, lb. or kg. lockout, and more

Lightweight and Easy to Carry: Built-in handle and recessed storage cavity securely holds the readout during storage and extended transport

Weight Measurements: In pounds or kilograms

PERFECT FOR WEIGHING
THE ELDERLY
&
HEAVY SET!



Befour Inc.
102 Progress Drive
Saukville, WI 53080

Phone
800-367-7109

Web
www.befour.com

Fax
262-284-5966

E-mail
sales@befour.com