

# BEFOUR

## PRECISION SCALES



Proudly Made in the USA

### PORTABLE BATTERY

Quality, ease of use, and versatility make this an all round value.

## Model PS-5701 Pro-BMI Scale

### Hi-Capacity PRO-BMI scales designed with heavy-weighing in



Battery power—no outlet or re-charging needed!!

The PS-5701 is truly a versatile scale that is ideal as a weigh-in scale where precise weight measurements are required. Its built-in body mass index calculator, makes it equally great for home health care professionals, athletic trainers, and weight loss clinics.

Convenient alkaline battery power allows user to take the scale anywhere - no outlet or re-charging needed!



#### MODEL PS-5701

Platform Dimensions .....	13" x 13" x 2"
Weight Capacity .....	500 lbs. or 225 kg.
Accuracy / Resolution .....	± 0.1 lb. (± 0.1 kg.)
Power Source .....	6 C-cell batteries
Readout .....	1" LCD digits
System Weight .....	12 lbs.
Warranty .....	3 Years

### FEATURES/SPECIFICATIONS:

**Pinpoint Accuracy:** ± 0.1 lb

**High Capacity:** 500 lb

**Programmable Settings:** include auto-off time, "Always-on" option, audible beep, lb. or kg. lockout, and more.

**BMI Calculator:** Operates in both user and assistant modes.

**Auto-zero, weight recall and reweigh**

**Plug & Weigh:** Exclusive cable system can be easily replaced if it is damaged.

**MotionLock Technology:** Accurate weight displays even with unsteady users.

**Push Button Calibration:** Set calibration via keypad entry with no tools needed.

### DISPLAY FEATURES

#### TARE CONSOLE



- Alkaline batteries provide over 80,000 weighings
- BMI calculator with programmable clothing allowance
- "Always-On" option
- Weight reads in pounds or kilograms
- Programmable pound/kilogram lockout
- Programmable auto-off: 12-seconds to 25-minutes
- Auto zero
- Weight recall and reweigh
- Audible sound indicates weight is "locked-in"
- Durable anodized aluminum construction
- Large 1" easy-to read digits
- 3-segment battery life indicator

### Model SC-1816 Soft-Sided Carry Case



Padded nylon soft case. Safely stores scale platform, readout, console, and cables. Great for storage, protection, and transport

### Model HC-1825 Hard-Shell Carry Case



Heavy-duty polyethylene hard-shell case with custom-cut foam interior is designed for added safety during transport

### OPTIONAL CARRY CASES:



Befour Inc.

102 Progress Drive Phone Web Fax E-mail  
Saukville, WI 53080 800-367-7109 www.befour.com 262-284-5966 sales@befour.com