


























BEFOUR

PRECISION SCALES

SCALE COMPARISON GUIDE

...the important information you need at a glance

	 PS-5700 500 pounds	 PS-6700 500 pounds	 PS-6600 ST 500 pounds	 PS-7700 900 pounds	 VS-0800 600 pounds
Weight Capacity	500 pounds	500 pounds	500 pounds	900 pounds	600 pounds
Platform Size (width x depth)	13" x 13"	16" x 18"	16" x 18"	20" x 24"	22" x 42"
Scale/System Weight	12 pounds	15 pounds	16 pounds	25 pounds	46 pounds
Accuracy/Resolution	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound
Display	1" LCD Digits	1" LCD Digits	3/4" LED Digits	1" LCD Digits	1" LCD Digits
Style	Portable	Portable	Portable	Portable	Portable
Warranty	3 Years	3 Years	3 Years	3 Years	3 Years
MotionLock					
Alkaline Battery Powered	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Rechargeable Batteries			<input type="checkbox"/>		
AC Plug In	Optional	Optional	<input type="checkbox"/>	Optional	Optional
Body Mass Index (BMI)				<input type="checkbox"/>	
Plug & Weigh					
	 FS-0900 500 pounds	 FS-0961 500 pounds	 PS-8070 500 pounds	 WH-1061 650 pounds	 WH-1070 650 pounds
Weight Capacity	500 pounds	500 pounds	500 pounds	650 pounds	650 pounds
Platform Size (width x depth)	16" x 18"	16" x 18"	16" x 18"	16" x 18"	16" x 18"
Scale/System Weight	17 pounds	17 pounds	25 pounds	17 pounds	25 pounds
Accuracy/Resolution	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound
Display	3/4" LED Digits	1" LCD Digits	1" LCD Digits	1" LCD Digits	1" LCD Digits
Style	Stand-on	Stand-on	Tilt & Roll Handrail	Stand-on	Tilt & Roll Handrail
Warranty	3 Years	3 Years	3 Years	3 Years	3 Years
MotionLock					
Alkaline Battery Powered		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Integrated Height Rod				<input type="checkbox"/>	<input type="checkbox"/>
AC Plug In	<input type="checkbox"/>	Optional	Optional	Optional	Optional
Body Mass Index (BMI)		<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>

Befour Inc.
102 Progress Drive
Saukville, WI 53080

Phone
800-367-7109

Web
www.befour.com

Fax
262-284-5966

E-mail
sales@befour.com



Proudly Made in the USA