

## Scale Comparison Guide - PS-6615 / PS-6600ST / PS-6700 /PS-6701

PS-6615



PS-6600ST



PS-6700



PS-6701



**PS-6615** with LCD digits.  
Console has a bracket to sit upright on a table top or can be mounted to the wall



**PS-6600ST** with LED digits.  
Console is wedge shape to sit upright on a table top.



**PS-6700 / PS6701** with LCD digits



Tare Only



Body Mass Index (BMI)

	PS-6615	PS-6600ST	PS-6700	PS-6701
Weight capacity	500 lbs.	500 lbs.	500 lbs.	500 lbs.
Platform Size (W x D )	16" x 18"	16" x 18"	16" x 18"	16" x 18"
Scale/System weight	15 lbs.	16 lbs.	15 lbs.	15 lbs.
Accuracy/Resolution	± 0.1 lbs.	± 0.1 lbs.	± 0.1 lbs.	± 0.1 lbs.
Display Console	1" LCD digits	1" LED digits	1" LCD digits	1" LCD digits
Power Source	4 D-Cell Alkaline Batteries	Rechargeable batteries or 120V AC	6 C-Cell Alkaline Batteries	6 C-Cell Alkaline Batteries
Warranty	3 years	3 years	3 years	3 years
Weighings	Over 100,000	1000+weighings per charge	80,000	80,000
MotionLock	<b>MOTIONLOCK</b>		<b>MOTIONLOCK</b>	<b>MOTIONLOCK</b>
Alkaline battery power	•		•	•
Body Mass Index (BMI)	•			•
Optional Module/App	•			

**OPTIONAL:** PS-6615 only, Bluetooth Smart Plug-in Module with Wrestling Weight Manager App

**BTM 415 Module**

**App**



**BEFOUR**  
WRESTLING  
WEIGHT MANAGER

**\*\*Please note:** If ordering the PS-6615, the Bluetooth Smart plug-in module is optional and can be ordered separately at a later date.