

Scale Comparison Guide - PS-6615 / PS-6600ST

PS-6615

PS-6600ST



PS-6615 with LCD digits.
Console has a bracket to sit upright on a table top or can be mounted to the wall



	PS-6615	PS-6600ST
Weight capacity	500 lbs.	500 lbs.
Platform Size (W x D)	16" x 18"	16" x 18"
Scale/System weight	15 lbs.	16 lbs.
Accuracy/Resolution	± 0.1 lbs.	± 0.1 lbs.
Display Console	1" LCD digits	1" LED digits
Power Source	4 D-Cell Alkaline Batteries	Rechargeable batteries or 120V AC
Warranty	3 years	3 years
Weighings	Over 100,000 + weighings	80,000+weighings per charge
<i>MOTIONLOCK</i>	<i>MOTIONLOCK</i>	
Alkaline battery power	•	
Body Mass Index (BMI)	•	
Optional Module/App	•	

PS-6600ST with LED digits.
Console is wedge shape to sit upright on a table top.



OPTIONAL: PS-6615 only, Bluetooth Smart Plug-in Module with Wrestling Weight Manager App

BTM 415 Module

APP



BEFOUR
WRESTLING
WEIGHT MANAGER

****Please note:** If ordering the PS-6615, the Bluetooth Smart plug-in module is optional and can be ordered separately at a later date.