

BEFOUR

SPORTS & FITNESS SCALES

Digital Scale with
Integrated Height Rod

Model WH-1061

Befour's WH-1061 weight and height measurement scale features a full-span, integrated digital height rod that quickly determines weight, height and BMI. With a weight capacity of 650 lbs., +/- 0.1lb accuracy, and MotionLock technology, which provides an accurate weight of unsteady users, this scale is an exceptional value.

The incredible battery life provides years of use on a single set of batteries—no outlet or recharging needed.

FEATURES & BENEFITS:

Integrated digital Height rod:

Simply raise the rotating headpiece and height automatically displays with BMI

Pin point accuracy of ± 0.1 lb.

High capacity: 650 lbs.

Body Mass Index:

Automatically calculates and displays BMI

Push Button Calibration:

Set calibration via keypad entry without the need of tools

Rugged and durable all-metal construction



Accurately calculates the weight of an unsteady user



SPECIFICATIONS:

Height Rod Range.....	38-88" (97-225cm)
Platform Dimensions.....	16" x 18" x 1 3/4"
System Weight.....	17 lbs
Resolutions.....	0.1 lb (0.05 kg)
Accuracy.....	± 0.1 lb (0.05 kg)
Power Source.....	4 D-cell Alkaline Batteries
Display Readout.....	1" LCD digits
Warranty.....	3 years

DISPLAY CONSOLE FEATURES



- Displays weight, height and BMI on a single screen
- **Programmable User Settings:** can be set to auto-off (12 seconds to 25 minutes), always-on, etc.
- **LB/KG Lockout:** Change setup to display only your desired unit of measure
- 3-segment battery life indicator
- Large 1" LCD digits
- No AC outlet or recharging needed (batteries last all season)
- Auto zero, recall, and reweigh



Height rod headpiece rotates 360° so the headpiece can be positioned away from patient's face when being unfolded and then raised above the patient's head.



Befour, Inc.

102 Progress Drive · Saukville, WI 53080 · P: 800-367-7109 · F: 262-284-5966
Befour.com · sales@befour.com