



















Weight Capacity	PS-5700 500 pounds	PS-6700 500 pounds	PS-6600ST 500 pounds	PS-7700 900 pounds	VS-0800 600 pounds
Platform Size (width x depth)	13" x 13"	16" x 18"	16" x 18"	20" x 24"	22" x 42"
Scale/System Weight	12 pounds	15 pounds	16 pounds	25 pounds	46 pounds
Accuracy/Resolution	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound
Display	1" LCD Digits	1" LCD Digits	0.4" LED Digits	1" LCD Digits	1" LCD Digits
Style	Portable	Portable	Portable	Portable	Portable
Warranty	3 Years	3 Years	3 Years	3 Years	3 Years
MotionLock					
Alkaline Battery Powered	◆	◆		◆	◆
Rechargeable Batteries			◆		
AC Plug-In	Optional	Optional	◆	Optional	Optional
Body Mass Index (BMI)				◆	
Connection (platform to display)					

					
Weight Capacity	FS-0900 500 pounds	FS-0961 500 pounds	PS-8070 500 pounds	WH-1061 650 pounds	WH-1070 650 pounds
Platform Size (width x depth)	16" x 18"	16" x 18"	16" x 18"	16" x 18"	16" x 18"
Scale/System Weight	17 pounds	17 pounds	25 pounds	17 pounds	25 pounds
Accuracy/Resolution	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound	± 0.1 pound
Display	0.4" LED Digits	1" LCD Digits	1" LCD Digits	1" LCD Digits	1" LCD Digits
Style	Stand-on	Stand-on	Tilt & Roll Handrail	Stand-on	Tilt & Roll Handrail
Warranty	3 Years	3 Years	3 Years	3 Years	3 Years
MotionLock					
Alkaline Battery Powered		◆	◆	◆	◆
Integrated Height Rod				◆	◆
AC Plug-In	◆	Optional	Optional	Optional	Optional
Body Mass Index (BMI)		◆		◆	◆

Befour, Inc.

102 Progress Drive • Saukville, WI 53080
800-367-7109 • Fax: 262-284-5966
befour.com • sales@befour.com



Scales proudly
made in the USA