

# BEFOUR

PRECISION SCALES



Proudly Made in the USA


## ANIMAL SCALE

Quality, reliability, and easy-to-use controls.

### Model VS-0800 Animal Scale

High capacity and extremely accurate...

BATTERY POWER -  
NO OUTLET  
OR RECHARGING  
NEEDED!



Our VS-0800 vet scale provides the most flexibility ever! Powered by ( 4 ) D-cell alkaline batteries, this scale can be used anywhere... on the floor, table, or as a mobile unit.

- Stainless steel cover and slip resistant mat are removable for easy cleaning
- 600 pound capacity provides all your weighing needs
- Durable all-metal construction
- Pin point accuracy of  $\pm 0.1$  lb.
- Console mounts on a wall, angles on a counter, or lays flat (brackets included)
- Exclusive "Plug & Weigh" cable system connects the platform to the console. Should the cable ever get damaged, replacement takes only minutes and can be easily done by the scale user.
- Industry-best 3 year warranty

### FEATURES/SPECIFICATIONS:

**MotionLock Technology:** Get an incredibly accurate measurement even with a moving animal

**Dual Tare Operation:** Automatic and keypad tare features help determine net weight of an animal

**Stored Tare Operation:** Easily deduct weight of a standard cage

**Push-Button Calibration:** Set calibration via keypad entry without the need of tools

**Industry-best 3 year warranty**

#### MODEL VS-0800

Platform Dimensions.....	42" x 22" x 1 $\frac{3}{4}$ "
Weight Capacity.....	600 lbs. or 270 kg.
Accuracy / Resolution .....	$\pm 0.1$ lb. ( $\pm 0.1$ kg.)
Power Source.....	6 C batteries (included)
Readout.....	1" LCD digits
System Weight .....	46 lbs.
Warranty.....	3 Years

### LCD DISPLAY FEATURES:



- Alkaline batteries provide over 80,000 weighings
- "Always-On" option
- Weight reads in pounds or kilograms
- Programmable pound/kilogram lockout
- Programmable auto-off: 12 seconds to 25 minutes
- Auto zero
- Weight recall and reweigh
- Audible sound indicates weight is "locked-in"
- Durable anodized aluminum construction
- Large 1" easy-to read digits
- 3-segment battery life indicator



Befour Inc.

102 Progress Drive  
Saukville, WI 53080

Phone  
800-367-7109

Web  
www.befour.com

Fax  
262-284-5966

E-mail  
sales@befour.com