



Proudly Made in the USA

WEIGHT + HEIGHT SCALE

Rugged construction, effortless portability...

Model WH-1070

Tilt and Roll Digital Scale with Integrated Height Rod!



WH-1070 is a truly versatile scale that is ideal for most weighing needs. Convenient battery power allows a user to place the scale anywhere - no outlet or re-charging needed! Sturdy, integrated handrail and digital height rod are excellent features providing stability for unsteady users and an accurate height and weight measurement.

Transport is easy ... simply tilt the scale back on its heavy-duty wheels and roll to the next location.

Height rod headpiece rotates 360° so the headpiece can be positioned away from the patient's face when being unfolded and then raised above the patient's head.



MOTIONLOCK

FEATURES/SPECIFICATIONS

- Huge 650 pound capacity with ± 0.1 lb. accuracy.
- Integrated digital height rod - simply raise the headpiece and height automatically displays with BMI.
- Ease of tilt and roll mobility.
- MotionLock Technology provides accurate measurement even for unsteady users.
- Standard D-cell alkaline batteries provide up to 100,000 weighings.
- Large platform for greater stability.
- Programmable to your specific needs.
- Rugged, lightweight, all-metal construction for years of service.
- Easy assembly.

MODEL WH-1070

Height Rod Range	38-88"/97-225cm
Platform Dimensions	16" x 18" x 1 3/4"
Weight Capacity	650 lbs. or 295 kg.
Accuracy/Resolution	± 0.1 lb. (± 0.1 kg.)
Power Source	4 D-cell batteries
Readout.....	1" LCD digits
System Weight.....	25 lbs.
Warranty.....	3 Years

WH-1070 DISPLAY FEATURES



- Easy keypad operation
- Unique display will show only the information you want
- Digital height measurement, simply press HEIGHT and raise the foldable headpiece.
- 3 segment battery-life indicator
- Programmable user settings can tailor console to your needs
- Push-button calibration
- Tare button is great for net weight

