



Proudly Made in the USA

Quick, reliable and extremely accurate
in high volume environments...

Model WH-1061 Digital Scale with Integrated Height Rod!



Befour's WH-1061 weight and measurement scale features a full span integrated digital height rod that quickly determines weight, height and BMI. With a weight capacity of 650 lbs., and +/- 0.1 lb. accuracy, the scale's MotionLock technology can weigh even the most unsteady users.

Height rod headpiece rotates 360° so the headpiece can be positioned away from the patient's face when being unfolded and then raised above the patient's head.



FEATURES/SPECIFICATIONS

- Huge 650 pound capacity with ± 0.1 lb. accuracy.
- Integrated digital height rod - simply raise the headpiece and height automatically displays with BMI.
- MotionLock Technology provides accurate measurement even for unsteady users.
- Standard alkaline batteries provide non-hazardous power.

Model WH-1061	
• Capacity:	650 lbs. (295 kg.)
• Platform Dimensions:	16"W x 18"D x 1 3/4"H
• Accuracy/Resolution:	± 0.1 lb. / (0.05 kg.)
• Height Measurement Range:	38-88" / 97-224 cm
• System Weight:	17 lbs.
• (4) D-cell batteries (included):	Provides up to 100,000 weighings
• Unique keypad:	Shows weight only, or weight, height, and BMI on a single screen
• Warranty:	Industry leading 3 years
• Made in the USA	

WH-1061

DISPLAY FEATURES



- Easy keypad operation
- Unique display will show only the information you want
- Digital height measurement, simply press HEIGHT and raise the foldable headpiece.
- 3 segment battery-life indicator
- Programmable user settings can tailor console to your needs
- Push-button calibration
- Tare button is great for net weight

