

BEFOUR

PRECISION SCALES



Proudly Made in the USA

WEIGHT + HEIGHT SCALE

Quick, reliable and extremely accurate in high volume environments...

Model WH-1061 Digital Scale with Integrated Height Rod!



Befour's WH-1061 weight and measurement scale features a full span integrated digital height rod that quickly determines weight, height and BMI. With a weight capacity of 650 lbs., and +/- 0.1 lb. accuracy, the scale's MotionLock technology can weigh even the most unsteady users.

Height rod headpiece rotates 360° so the headpiece can be positioned away from the patient's face when being unfolded and then raised above the patient's head.



FEATURES/SPECIFICATIONS

- Huge 650 pound capacity with ± 0.1 lb. accuracy.
- Integrated digital height rod - simply raise the headpiece and height automatically displays with BMI.
- MotionLock Technology provides accurate measurement even for unsteady users.
- Standard alkaline batteries provide non-hazardous power.

MODEL WH-1061

Height Rod Range.....	38-88"/97-225cm
Platform Dimensions	16" x 18" x 1 3/4"
Weight Capacity	650 lbs. or 295 kg.
Accuracy/Resolution	± 0.1 lb. (± 0.1 kg.)
Power Source.....	4 D-cell batteries
Readout.....	1" LCD digits
System Weight	17 lbs.
Warranty.....	3 Years

WH-1061

DISPLAY FEATURES



- Easy keypad operation
- Unique display will show only the information you want
- Digital height measurement, simply press HEIGHT and raise the foldable headpiece.
- 3 segment battery-life indicator
- Programmable user settings can tailor console to your needs
- Push-button calibration
- Tare button is great for net weight



Befour Inc.
102 Progress Drive
Saukville, WI 53080

Phone
800-367-7109

Web
www.befour.com

Fax
262-284-5966

E-mail
sales@befour.com